




Dragon – Establishing ECMC Profile

To establish an ECMC Dragon Profile only takes a few minutes. When you first log into Dragon you will be prompted to create the profile by reading a few short scripts. The screens are easy to follow, but the steps are outlined below. If you have any issues, please contact the Help Desk at 898-4477.

Microphone Tips

- Point the microphone element to the side of or below your mouth to avoid noise from breathing. Do not hold the microphone directly in front of your mouth.
- Hold the microphone in a comfortable position.
- Make sure that the microphone is about 1 to 3 inches away from your mouth. You should keep the microphone distance in this range, but you don't have to always keep the microphone at exactly the same distance.
- Hold the microphone so that the front points toward your mouth. Don't talk into the side of the microphone.



Step by Step to Create Profile

1.	 <p>Open Dragon</p>		<p>If you need to log in, use your ECMC Network user name and password</p>
2.			<p>Click Next></p>

<p>3.</p>		<p>Select your age group to help Dragon customize your profile.</p>	<p>Click on Next></p>
<p>4.</p>		<p>Select:</p> <ul style="list-style-type: none"> • Region: United States • Vocabulary: pick your specialty or select General Medical 	<p>Click on Next></p>
<p>5.</p>		<p>If you have accented English, you may select an accent. If not, or if your accent is not listed, click on Standard.</p>	<p>Click on Next></p>

Dragon – Establishing ECMC Profile

<p>6.</p>		<p>The microphone that is plugged in will be defaulted here.</p> <p><i>ECMC is providing PowerMic II.</i></p>	<p>Click on Next></p>
<p>7.</p>		<p>Review the selections you made.</p> <p>If you need to correct anything, click on <Back.</p>	<p>If everything is correct click on <input checked="" type="checkbox"/> I have read and agree... and click on Create</p>
<p>8.</p>	<p>You will receive a message that says, “Your profile is being created”</p>	<p>This may take a few minutes.</p>	
<p>9.</p>		<p>You will read a few short paragraphs so that Dragon can begin learning your speech patterns.</p>	<p>When ready, click on Next></p>

<p>10.</p>		<p>Position your microphone correctly <i>Tips are listed at the top of this document</i></p>	<p>Click on Start Volume Check and begin reading the paragraph on the screen.</p>
<p>11.</p>	<p>When you have successfully read the paragraph, you will hear a beep and Volume has been adjusted successfully will appear at the top of the screen.</p>		<p>Click on Next></p>
<p>12.</p>		<p>You will need to read one more paragraph</p>	<p>Position the mic and click on Start Quality Check.</p>
<p>13.</p>	<p>Read the paragraph. When Dragon has enough information, it will beep and it will display Check Microphone is complete at the top of the screen.</p>		<p>Click on Next></p>
<p>14.</p>	<p>Your profile is being updated will appear</p>	<p>This step may take a few minutes. When done, it will automatically advance you and turn on Dragon.</p>	

Dragon – Establishing ECMC Profile

<p>15.</p>		<p>You may click on the links for more information</p>	<p>Click on Finish</p>
<p>16.</p>	<p>At the top of your screen you will see the Dragon toolbar.</p>		
<p>17.</p>	<ul style="list-style-type: none"> • Open up the program that you would like to dictate into. • Press and hold the red dictate button on the mic to record your speech 		